



Mass Spectrometry for Nanovectors Characterization



Aura TINTARU


Aix*Marseille
université







Outline

- **Introduction to Mass Spectrometry**
 - Ionization Sources
 - Analyzers
- **Tandem Mass Spectrometry**
- **Ion Mobility Mass Spectrometry**
 - Principles et Definitions
 - Instrumentation
 - Applications
- **Conclusions**





Mass Spectrometry

Distinguish between ionized molecules depending on...




... mass
m
(isotopes)





... charge
z
(ions)

m/z





Mass Spectrometry






Nano2Clinic
CA17140


Mass Spectrometry




MS spectrum



Electrospray Ionization (ESI)




- This is often coupled with LC (LC-MS)
- The solvated analyte is fed through a capillary tube with very high charge differential
- The charged droplet evaporates and ionizes the analyte




High Resolution Mass Spectrometry

- ◆ Fast method to determine the elementary composition (error <5ppm)
- ◆ Use of low quantities
- ◆ Allows analysis of all types of organic molecules





Nano2Clinic
CA17140

High Resolution Mass Spectrometry




Theoretic: m/z 587.1019





Nano2Clinic
CA17140

High Resolution Mass Spectrometry



Intensity,



2000




Theoretic: m/z 587.1019

Experimental: m/z 587.1011


Error: -1.4 ppm






Nano2Clinic
CA17140


Tandem Mass Spectrometry






Tandem Mass Spectrometry

- ◆ Useful tool for structural characterization
- ◆ Appropriate for structural characterization of highly symmetrical molecules: polymers, dendrimers, sugars
- ◆ Suitable for structural anomalies screening



Nano2Clinic
CA17140


Dendrimers as Drug Delivery Systems




- Highly symmetrical structures
- Intrinsic cavities
- Functionalized surface



Vectorization agents



Dendrimers as Drug Delivery Systems




guest

: DNA, RNA, drugs, ...

- Highly symmetrical structures
- Intrinsic cavities
- Functionalized surface



Vectorization agents



Nano2Clinic
CA17140


Dendrimers as Drug Delivery Systems

DENDRIMERS


Complex molecules




Original analytical approaches
adapted to each specific topic




Main goal: Accurate characterization




Dendrimers as Drug Delivery Systems




Goal: Detection and characterization of the defective molecules





Nano2Clinic
CA17140


MS Study of nanovectors



→ Complementary approach to non-covalently dendrimer-drug systems

Drug conjugation on the dendrimer arms


MS Study of nanovectors



~~Triazole cycle = linker~~


Should play a role in the drug delivery process

Drug



Triazole nucleoside = anticancer activity

- High drug loading
 - Can be much easier to control and predict in complex media
 - Overcome the problem of poor solubility of triazole nucleoside
 - Preserve the tumor targeting properties of nanosystems



Nano2Clinic
CA17140

MS Study of nanovectors

Intensity (a.u.)

5.0e4

3.0e4

2.0e4

1.0e4

500

700

900


1100


1300

1500

1700


m/z






Nano2Clinic
CA17140


MS/MS Study of nanovectors






Nano2Clinic
CA17140


MS/MS Study of nanovectors






Nano2Clinic
CA17140


MS/MS Study of nanovectors






Nano2Clinic
CA17140


MS/MS Study of impurity A






Nano2Clinic
CA17140


MS/MS Study of impurity A








Nano2Clinic
CA17140


MS/MS Study of impurity A






MS Study of nanovectors







Nano2Clinic
CA17140

MS/MS Study of nanovectors






MS/MS Study of nanovectors




Fragmentation rules → used for the structural characterization of the two impurities B and C




Nano2Clinic
CA17140

MS/MS Study of nanovectors


Impurity B



Impurity C




Structure of the “defective” molecules




Nano2Clinic
CA17140

MS/MS Study of nanovectors





Amphiphilic Dendrimers



Main goal: Highest sample purity required for any for biological application

Nanotheranostics






Imaging

PET

Positron Emission Tomography


Gamma rays emitted indirectly by a positron-emitting radioligand



SPECT

Single Photon Emission Computed Tomography


Gamma-emitting radioisotope




MRI

Magnetic Resonance Imaging

Magnetic relaxation of water and fat tissue





Imaging

PET

Positron Emission Tomography


Gamma rays emitted indirectly by a positron-emitting radioligand



SPECT

Single Photon Emission Computed Tomography

Gamma-emitting radioisotope




MRI

Magnetic Resonance Imaging


Magnetic relaxation of water and fat tissues




Amphiphilic Dendrimers



Versatile Goal: Accurate characterization of EACH KIND of micelle







Nano2Clinic
CA17140

NOTA-Amphiphilic Dendrimers

MS & MS/MS Study

Expected: $C_{111}H_{222}N_{28}O_{34}$








NOTA-Amphiphilic Dendrimers

MS & MS/MS Study

Expected: $C_{111}H_{222}N_{28}O_{34}$






Nano2Clinic
CA17140


NOTA-Amphiphilic Dendrimers

MS & MS/MS Study

Expected: $C_{111}H_{222}N_{28}O_{34}$


Detected: $[C_{111}H_{222}N_{28}O_{34} + 3H]^{3+}$






Nano2Clinic
CA17140

Ga^{III}-NOTAAmphilic Dendrimers





39




Nano2Clinic
CA17140

Ga^{III}-NOTAAmphilic Dendrimers



Expected: $C_{111}H_{222}N_{28}O_{34} Ga^{(III)}_4$


40




Nano2Clinic
CA17140

Ga^{III} -NOTAAmphilic Dendrimers

PET experiments









Migration of the dendrimer in the cancerous tissue after 10min



Ion Mobility: a New Dimension for MS

IMS separates ions based on their mobility through a neutral gas

The **mobility** of an ionized molecule depends on ...




... its size

... its charge

... its shape


already discriminated in MS




Ion Mobility: Principle




Which will "land" first?



Ion Mobility: Principle









Ion Mobility: Principle

In IMS, ions are separated as a function of their mobility through an inert gas under the effect of an electric field.


Same masses
but different
conformations



- Mobility of an ion depends on its WEIGHT, CHARGE STATE & SHAPE




Ion Mobility: Principle




Drift times (ms)


Ion Size (\AA^2)



Collision Cross Section (CCS)





IMS: Gas phase diffusion




Nano2Clinic
CA17140


Ion Mobility: Instrumentation







Ion Mobility: Instrumentation






Nano2Clinic
CA17140


Ion Mobility: Instrumentation






Nano2Clinic
CA17140


Ion Mobility: Instrumentation






Nano2Clinic
CA17140



Ion Mobility: Instrumentation






Nano2Clinic
CA17140

Ion Mobility – Mass Spectrometry




Triple correlation: m/z , drift time and abundance
Collision Cross Section (CCS)






Ion Mobility – Mass Spectrometry

- ◆ Distinguish between different conformations of the same molecule
- ◆ Isomers separation
- ◆ Discrimination between isobaric species with different shape




Isomers Separation





$\Delta\Omega < 3\text{\AA}^2$



Isobaric amino acids differentiation




Isomers Separation




Isomers Separation




IMS Applications Dinucleotides




dTA






Nano2Clinic
CA17140


IMS Applications Peptides: Bradikinin



Water: Methanol






Water: Dioxane



Nano2Clinic
CA17140


IMS Applications Proteins: Ubiquitin






Nano2Clinic
CA17140


IMS – MS Study of Dendrimers



PAMAM
Globular Shape



PAMAM
Fan Shape




Nano2Clinic
CA17140


IMS – MS Study of Dendrimers

PAMAM Fan Shape


1+



2+



3+



Theoretical
Experimental

221 \AA^2


219 \AA^2

204 \AA^2

208 \AA^2



241 \AA^2


229 \AA^2



Nano2Clinic
CA17140

In^{III}-NOTA Amphiphilic Dendrimers







Nano2Clinic
CA17140

In^{III}-NOTA Amphiphilic Dendrimers


IMS-MS Study



IM-MS → CCS calculation (\AA^2)



	NOTA	Dendrimer-NOTA
free	149.3	581.2
Ga (III)	154.5	580.7
In (III)	158.1	570.9

NOTA "cage" folded




Nano2Clinic
CA17140

In^{III}-DOTA Amphiphilic Dendrimers




SPECT experiment

64




IMS – MS/MS Study of Dendrimers




- Micelles-like compounds → highly medium dependent
- High flexibility
- Larger intrinsic cavities

} novel drug delivery systems




Nano2Clinic
CA17140

IMS – MS/MS Study of Dendrimers




Dendrimer or polymer MS/MS behavior ?




Nano2Clinic
CA17140


IMS – MS/MS Study of Dendrimers




Polymer Fragmentation




Dendrimer Fragmentation






IMS – MS/MS Study of Dendrimers






Polymer Fragmentation



Dendrimer Fragmentation


IMS – MS/MS Study of Dendrimers



Experimental



CCS : 441 Å²

CCS : 498 Å²



Nano2Clinic
CA17140


IMS – MS/MS Study of Dendrimers




Theoretical / Experimental

CCS : 445 Å² / 441 Å²

CCS : 515 Å² / 498 Å²






Conclusions

Mass Spectrometry powerful technique:

- ◆ **Accurate Mass Measurements**
- ◆ **Structural Characterization**
- ◆ **Isomers Separation**
- ◆ **Differentiation Based on Ions Shape**



Nano2Clinic
CA17140

Acknowledgments



Dr. Ling PENG

Dr. Gilles QUELEVER

Dr. Zhenbin LYU

Dr. Mei CONG

Dr. Jingjie TANG

Dr. Adela HUANG



Dr. Christophe CHENDO

Dr. Valérie MONNIER

Dr. Fabio ZIARELLI



UNIVERSITÀ
DEGLI STUDI DI TRIESTE

Pr. Sabrina PRICL

Pr. Erik LAURINI



